VIRUS DISINFECTION

Durr Dental's hygiene range destroys all viruses, both enveloped and nonenveloped, such as polio and the Norovirus. Included in the range are solutions for disinfection of instruments as well as surfaces. Its ID213 Instrument Disinfection can impressively protect instruments for up to a fortnight. It is also highly compatible by virtue of its corrosion inhibitors.

For cleaning small areas there's FD333 Disinfection Wipes, a scent and residue free, fast drying wipe that is fully virucidal.

For larger areas there's FD333 Surface Disinfection. This high yielding solution works fast and the formulation is safe to use on a range of materials.

For more information call 01536 526740.



Reader Enquiry – 013

EFFICACY AND SAFETY

The DUX Dental team will be on hand at The Dentistry Show to demonstrate the array of simple to use and effective products. Image Alginate is a dust- and bubble-free alginate, designed for optimum absorption and maximum detail for impressions.

ZONEfree is a translucent and eugenol-free zinc oxide provisional cement.

PeelVue+ is a sterilisation pouche that ensure decontamination of instruments.

Bib-Eze is a lightweight, elasticised and disposable bib holder that protect patients from the potential infection risk.

To find out more and for free

samples visit DUX Dental on stand G92.

For more information email info@dux-dental.com or visit www.dux-dental.com









Reader Enquiry – 014

SAVING TEETH

EndoCare aims to become a valuable member of the team, helping and assisting with aspects of endodontic treatment. It welcomes contact from referring practices, and the team is happy to discuss cases and offer support.

Endocare takes referrals for anything from post removals to advanced microsurgery techniques. It will treat anxious patients and those dealing with anaesthetic problems and other

issues.
When patients visit it aims

challenging

to treat them with warmth, compassion and honest empathy. A member of the team will always be on hand to answer any questions or address any concerns the patient may have.

For more information call 020 7224 0999 or visit www. endocare.co.uk



Reader Enquiry – 015

NEW AUTOCLAVE

Echmann's new autoclave, the LittleSister SES 3000B, has recently been trialled in 12 dental practices and facilities across the UK.

Decontamination technician at The University of Portsmouth Dental Academy, Michael Tyler, says: "This was great for the

students to have access to – the display interface was clear, it was user-friendly, with regards to maintenance, identifying potential errors was easy, and it

handled a full load effectively.

"The team at Eschmann are receptive and interested in feedback. We are looking forward to using the new autoclave with students."

For more information please

For more information please call 01903 753322 or visit www.eschmann.co.uk







Reader Enquiry - 016

MAXIMISE TAX CLAIMS

When filing tax returns, maximise tax deductions and claim against everything allowable to get the best possible tax savings.

When it comes to maximising tax claims and saving money, figurit can help. This web based accountancy and tax solution for dental associates can take care of the accurate and timely filing of tax returns and bring the most in

tax savings.

With specialist
knowledge, it can
provide a complete
list of dental specific
claimable expenses and ensure

information is provided on how to record these accurately.

Because figurit knows both the dental industry and tax laws, everything claimed is wholly allowable.

For more information visit www.figurit.co.uk



Reader Enquiry – 017

REINFORCE RESTORATIONS

Claim a free sample of GC's everX Posterior by completing the online form on GC's Facebook page.

GC launched the 'adventures of everX' e-comic on Facebook to highlight the strengths of everX and to show how it can save a practice money.

GC everX Posterior is a fibre-reinforced composite designed to be used as dentine replacement. The fibres make it a sub-structure to reinforce composite restoration in large cavities. Fibres will prevent and stop crack propagation through the filling. It opens possibilities for restorations of extensive cavities at chairside and is the

answer to the growing demand for a restorative alternative. For more information call 01908 218 999.



Reader Enquiry – 018