MARKETPI ACE

DENTAL SKY

Dental Sky's CMA System was designed to simplify endodontic treatments by minimising the number of instruments necessary and providing one uncomplicated sequence for both treatment and retreatment.



CMA stands for coronal, median and apical and provides a simple sequence with only four core instruments in the range: coronal c to flare out the coronal portion of the canal, median m to enlarge the middle part of the canal and apical fine a1 and apical a2 to enlarge the apical portion of the canal

CMA provides safe, secure and reliable rotary instruments with several special design features including a shorter handle to improve ease of access to molars and a helical shape to aid removal of debris from the canal. CMA instruments also have a cross-section with three cutting angles for excellent cutting efficiency and feature a non-cutting tip to ensure the instruments respect the trajectory of the root canal.

CMA files are available in starter kits in 21mm, 25mm or 29mm and in refill packs of six instruments. To download the full CMA product brochure, please visit the website.

0800 294 4700 www.dentalsky.com

FILHOL DENTAL

For root post and core build-up, Filpost is designengineered to be easier, faster and safer to place, while preserving more healthy tooth structure and being stronger in use

Filpost is made of 99.8% pure titanium and it is biocompatible and compatible with all dental materials that won't corrode.


It can be customised to suit the canal without risk of fracture, enabling easy insertion of multiple posts into converging canals. No drilling is required so avoiding the risk of perforation. Its anatomical shape minimises dentine removal. Retention grooves along the post that are working together with retention grooves (formed within the canal surface by the special universal groover), create a passive interlock, which strengthens retention.

Filhol Dental has specialised in posts for more than 30 years 01386 841 864 info@filhol.com www.filhol.com



PERIPRODUCTS

Dental professionals can email their practice logo today to Periproducts and pick up a free, customised Ultradex oral healthcare kit at the BDIA Dental Showcase at stand O03.



between professional and home care, increase patient compliance, enhance treatment outcomes and encourage patients to become an active participant in their own oral health.

The oral health care packs also offer outstanding value for money and are a convenient way to increase practice profile and revenue.


Backed by scientific research, the clinically proven technology (stabilised chlorine dioxide) within Ultradex protects against plaque, restores natural whiteness and eliminates bad breath instantly for at least 12 hours.

Send the practice details and logo today to dental@periproducts.co.uk and pick up a free Ultradex oral healthcare kit customised with the practice details at the BDIA from stand

dental@periproducts.co.uk www.periproducts.co.uk

SIRONA DENTAL SYSTEMS

Precise, fast and safe - intraoral radiography is now going mobile. Since January this year, The Dental Compan Sirona has been supplying users of the intraoral



More importance is placed on the flexibility of treatment devices and systems in modern dental practices; this is why Sirona made the intraoral X-ray tube Heliodentplus mobile, to offer even more flexibility for dentists. As well as being able to attach the X-ray tube on the wall or ceiling, users can now also attach it to a mobile stand. This flexibility not only improves the daily workflow, but also makes X-rays more comfortable for patients. This assembly option is suitable for clinics and nursing homes so even if patients are immobile X-rays that are necessary for treatment can still be taken without having to move the patient to a practice. In large practices where several dentists work together or that have multiple treatment rooms, the mobile stand allows Heliodentplus to be quickly and easily moved from room to room.

0845 071 5040 info@sironadental.co.uk

ENDOCARE

Endocare brings together a talented group of endodontists who are ready and available to assist dentists and their patients with all aspects of this challenging procedure.

Endocare's light and airy practices are located across London and are equipped with the latest technology to ensure patients receive the highest standard of care.

From the moment the patient arrives, there will always be a member of Endocare's team who will be able to answer questions and alleviate any concerns the patient might have

EndoCare

For complete peace of mind, Endocare will always communicate with the referring dentist at each stage of their patient's time at the clinic and will return them to their dentist pain-free, confident and happy with the quality of treatment they have received

020 7224 0999 www.endocare.co.uk

GC UK

GC's Gradia Core and Fiber Post system adds another dimension to GC's growing minimal intervention programme by allowing dentists to complete an aesthetic core and post luting with minimum tooth preparation.

With the smallest possible post size the remaining tooth structure is maximised, resulting in less stress and

avoidance of weakening of the root, which may be a cause of root canal treatment failure. Conveniently dispensed by way of an auto-mixing tip, Gradia Core is dispensed directly into the canal and direct core build-up. Being thixotropic the material will not slump, which is useful when used in the upper jaw; and as Gradia Core is a complete luting and core build-up material in one, the core can be quickly and easily formed in just one appointment, saving the dentist's and patient's time.

This dual-cured material is highly filled (75%) providing strong support for the prosthesis. Shrinkage is minimal to avoid stress factors and, being radiopaque Gradia Core and Fiber Posts, they are easy to detect on radiographs. GC's Fiber Posts contain a speical light transmission element to accelerate the curing that's particularly useful in deep canals.

01908 218 999

SIRONA DENTAL SYSTEMS

Sicat, a specialist for 3D planning and therapeutic solutions, is now presenting the market release of Sicat Function. This is a

solution that precisely merges diagnostic patient information received from a DVT

(Galileos)/CT, a jaw motion tracker (Sicat JMT+) and optical impressions (Cerec, Sirona). With convincing results; Sirona says this is the first time that real patient-individual movement of the lower jar can be visualised true to the anatomy and within the 3D volume.

The demand is tremendous: more than six million people suffer from craniomandibular dysfunctions. Up to now, the diagnosis of such craniomandibular dysfunctions has been complex and time consuming while conventional approaches require a high degree of experience.

Sicat Function offers the possibility to meet these high requirements through the presentation of more precise and patient-individual information. The Sicat JMT+ records all degrees of freedom and movements of the lower jaw in a highly precise

0845 071 5040 info@sironadental.co.uk

VOCO

The trusted Rebilda post system is now joined by the extended Rebilda post system set II, which additionally contains five posts of 1.0mm diameter, one drill size 1.0mm and a second reamer. The new set also contains 20 Singledose of the dual-curing universal adhesive Futurabond U.

The new set, just like the existing Rebilda post system set, increases clinical safety as all the components have been chosen to match each other and are ready at hand as a complete set to be used for post-endodontic treatments.

A coronal build-up can be accomplished in just a few, stress-free steps using Rebilda DC as the luting and build-up composite, Futurabond U, the dual-

curing universal adhesive, Rebilda Post, the glass fibre-reinforced composite root post with exactly fitting pilot and root canal drills as well as Ceramic Bond – a high strength silanisation agent

info@voco.com www.voco.com



