

Endodontics

Private Dentistry presents a selection of endodontic products and equipment

COLTENE WHALEDENT

To enhance its already wide endodontic range, Coltene has launched the new Endo Cleaning Pack, which contains the products needed for irrigation and drying of the root canal during endodontic treatment. Included are Surgitip-endodontic aspirator tips for effective canal accessibility and aspiration, Canalpro NaOCl 3% for irrigation and debridement, Roeko canal brush to remove loosened tooth matter, paper points and cotton pellets, Canalpro slotted end tips and colour-coded syringes for fast and easy delivery of solutions. The total saving of buying the Endo Cleaning Pack is more than 50%.



- 01444 235486 exts 223/224
- www.coltene.com

ENDOCARE

Endocare excels at diagnosing and treating dental pain. Its team of trained endodontists is committed to relieving pain and saving teeth. In the last 10 years, the Endocare team has saved more than 10,000 teeth from extraction.

Endocare's philosophy is simple: it aims to provide high quality endodontic treatments in a warm and caring environment. The patients' needs are always Endocare's number one concern.

Those who have a difficult diagnosis or need help with any aspect of root canal treatment, give Endocare's team a call today. It welcomes contact from all referring practices and are ready and willing to assist in any way they can. The team excels at treating nervous or anxious patients and other challenging issues such as hot pulp, microsurgery and retreatments.

- 020 7224 0999
- www.endocare.co.uk



EVIDENT


Up to 82% of dentists are reported to be suffering from musculo-skeletal complaints such as back or neck pain; it's a serious disorder and can have a negative impact on a dentist's quality of life.

Providing dental loupes and lights to the UK market for more than 20 years, Evident's range of Examvision loupes has been designed to allow the ideal posture to be attained, so dentists can start to see the bigger picture while maintaining a better working posture.

Whatever the loupe requirement, simply choose from the following: Kepler Kompakt – the powerful yet compact loupe in 3.5x, 4.6x and 5.7x, HD – the stylish and lightweight loupe available in 2.3x to 3.8x or Essential – the cost-effective 2.5x loupe.

Evident also provides Examvision lights to help improve visibility of the working site with the Essential, Focus, XP and the Nano that are compatible with Essential and premium batteries.

- 0500 321 111
- www.evident.co.uk



GC UK

GC's Gradia Core and Fiber Post system adds another dimension to GC's growing minimal intervention programme; it allows dentists to complete an aesthetic core and post luting with minimum tooth preparation.

With the smallest possible post size, the remaining tooth structure is maximised, with the result of less stress and an avoidance of weakening of the root that may be a cause of root canal treatment failure.

Conveniently dispensed by way of an auto-mixing tip, Gradia Core is dispensed directly into the canal and direct core build-up. Being thixotropic, the material will not slump, which is particularly useful when used in the upper jaw; and as Gradia Core is a complete luting and core build-up material in one, the core can be quickly and easily formed in just one appointment, saving dentists and patients time.

- 01908 218 999



VELOPEX

The ergonomically designed seat of this saddle chair directs body automatically in the natural, health-promoting posture.

Once the user has achieved this, it is easy to work all day. Velopex says the two-part seat is the 'only way' of sitting that is beneficial to the user's health.

Contact Velopex to book a free, in-practice demo. Dental professionals will be able to use the chair for an extended period of time to experience the postural and health benefits offered when using Velopex's saddle.

- 0208 965 2913
- dentistry@velopex.com



VOCO

Rebilda Post is a component of a complete, optimally coordinated post build-up system. As with Rebilda DC, the proven core build-up material, the endodontic post consists of a dimethacrylate matrix so that a reliable compound structure is attained under the build-up of a stable mono-block.

The complete endodontic post set Rebilda Post system increases clinical safety, since all components are coordinated with each other and all components are at hand if and when an endodontic post procedure is pending.

Rebilda Post system is adapted to the individual work steps of post cementation and subsequent core build-up. The Rebilda Post system contains five posts each of diameters 1.2mm, 1.5mm and 2mm, one drill each of matching sizes and one reamer of diameter 0.7mm. The new Rebilda Post System – set II now provides as an extra five posts of 1mm diameter, one drill size 1mm and an additional reamer. The root post, particularly in size 1.0mm, makes for substance-preserving and safe treatment of teeth with a narrow root canal.

- info-uk@voco.com
- www.voco.com

